# TANITA®

Monitoring Your Health

# *Body Composition Guide for InnerScan®*

**Eng**

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# WHO CAN USE A BODY COMPOSITION MONITOR

---

This Body Composition Monitor is intended for adults aged 18-99 years. Children aged 5-17 years can use the monitor for weight and body fat percentage readings and Healthy Range Indicator only ; the other features are not applicable to children.

Your monitor is also equipped with an athlete mode for adults with athletic body types. Tanita defines an athlete as a person involved in intense physical activity of approximately 10 hours per week and has a resting heart rate of 60 beats per minute or less. Individuals who have been an athlete for a number of years but currently exercises less than 10 hours per week can also use athlete mode.

Pregnant women should use the weight function only.  
All other functions are not intended for pregnant women.

The Tanita InnerScan Body Composition Monitor is intended for home use only. This unit is not intended for professional use including hospitals, medical or fitness facilities; it is not equipped with the adequate specifications required for heavy usage. Using the unit in a professional location will invalidate your warranty.

The Tanita InnerScan Body Composition Monitor product provides readings for informational purposes only. This product is not intended to diagnose or treated any disease or abnormalities. Please consult with your physician about any questions or concerns related to your health.

## WHY MONITOR BODY COMPOSITION?

Body composition monitors are designed for healthier living by giving you an insight into key health indicators that will enable you to monitor the impact of changes to your lifestyle:

- See the impact of a change in diet on the body to make sure you are dieting in the healthy way
- Fine-tune your fitness program by monitoring progress of muscle mass and Basal Metabolic Rate
- Monitor the level of visceral fat, which has been linked as a possible risk factor for developing Type 2 diabetes and/or heart disease.
- Set a target for your physique and monitor your progress towards it.

## HOW DOES A BODY COMPOSITION MONITOR WORK?

Tanita Body Composition Monitors calculate your body composition using Bioelectrical Impedance Analysis (BIA).

Safe, low-level electrical signals are passed through the body via the patented Tanita foot pads on the monitor platform. It is easy for the signal to flow through fluids in the muscle and other body tissues but meets resistance as it passes through body fat , as it contains little fluid. This resistance is called impedance. The impedance readings are then entered into medically researched mathematical formulas to calculate your body composition.

## WHEN IS THE BEST TIME TO USE MY BODY COMPOSITION MONITOR?

Your body water levels naturally fluctuate throughout the day and night. Any significant changes in body water may affect your body composition readings; for example, the body tends to be dehydrated after a long night sleep so if you take a reading first thing in the morning your weight will be lower and your body fat percentage higher. Eating large meals, drinking alcohol, menstruation, illness, exercising, and bathing may also cause variations in your hydration levels.

To get the most reliable reading it is important to use your Body Composition Monitor at a consistent time of day under consistent conditions. We suggest taking a reading before your evening meal.

## FOR SEGMENTAL BODY COMPOSITION MONITOR

Tanita has developed an accurate segmental body composition analyzer which correlates with the gold standard, DEXA, Dual Energy X-ray Absorptiometry. Tanita can now analyze the distribution of body fat as well as muscle in one's body.

Using this Segmental Body Composition Scale/Monitor, distribution in right arm, left arm, right leg, and left leg of body fat % and muscle mass can be determined. This is especially useful for anyone who is monitoring the balance of left and right side of the body or trying to build a particular part of one's body. It has been said that men tend to carry body fat in the upper body while women have a tendency to carry body fat in lower body. In addition, as you age, there is a tendency to lose muscle and gain fat.

It is important to monitor change in the distribution of body fat and muscle.

# WHAT IS BODY FAT PERCENTAGE?

(Applicable age 5-99)


## Body fat percentage is the amount of body fat as a proportion of your body weight.

Reducing excess levels of body fat has shown to reduce the risk of certain conditions such as high blood pressure, heart disease, diabetes and cancer.

The chart below shows the healthy ranges for body fat.

**Body Fat Ranges** for Standard Children <sup>1,2</sup>

**Body Fat Ranges** for Standard Adults <sup>3,4</sup>

- <sup>1</sup> *Body Fat Reference Curves for children Targeted at BMJ (British Medical Journal) Draft 1-AMP 19 June 2004 (by Dr Andrew)*
- <sup>2</sup> *Gallagher D et al. Am J Clin Nutr 2000, 72:694-701 Healthy percentage body fat ranges: an approach for developing guidelines based on body mass index.*
- <sup>3</sup> *Based on NIH/WHO BMI Guidelines.*
- <sup>4</sup> *As reported by Gallagher, et al, at NY Obesity Research Center. To determine the percentage of body fat that is appropriate for your body, consult your physician.*

## Healthy Range Indicator

Your Body Composition Monitor automatically compares your body fat percentage reading to the Healthy Body Fat Range chart. After your body fat percentage has been calculated, a black bar will flash along the bottom of the display, identifying where you fall within the Body Fat Ranges for your age and gender.

- (-)** : Underfat; below the healthy body fat range. Increased risk for health problems.
- (0)** : Healthy; within the healthy body fat percentage range for your age/gender.
- (+)** : Overfat; above the healthy range. Increased risk for health problems.
- (++)** : Obese; high above the healthy body fat range. Greatly increased risk of obesity-related health problems.

## WHAT IS TOTAL BODY WATER PERCENTAGE?

(Applicable age 18-99)

**Total Body Water Percentage is the total amount of fluid in a person's body expressed as a percentage of their total weight.**

Water plays a vital role in many of the body's processes and is found in every cell, tissue and organ. Maintaining a healthy total body water percentage will ensure the body functions efficiently and will reduce the risk of developing associated health problems.

Your body water levels naturally fluctuate throughout the day and night. Your body tends to be dehydrated after a long night sleep and there are differences in fluid distribution between day and night. Eating large meals, drinking alcohol, menstruation, illness, exercising, and bathing may cause variations in your hydration levels.

Your body water percentage reading should act as a guide and should not be used to specifically determine your absolute recommended total body water percentage. It is important to look for long-term changes in total body water percentage and maintain a consistent, healthy total body water percentage.

Drinking a large quantity of water in one sitting will not instantly change your water level. In fact, it will increase your body fat reading due to the additional weight gain. Please monitor all readings over time to track the relative change.

Every individual varies but as a guide the average total body water percentage ranges for a healthy adult are:

**Female** : 45 to 60%

**Male** : 50 to 65%

**Source** : Based on Tanita's Internal Research

**Note:** The total body water percentage will tend to decrease as the percentage of body fat increases. A person with a high percentage of body fat may fall below the average body water percentage. As you lose body fat, the total body water percentage should gradually move towards the typical range given above.

## WHAT IS VISCERAL FAT RATING?

(Applicable age 18-99)

**This feature indicates the rating of visceral fat in your body.**

Visceral fat is the fat that is in the internal abdominal cavity, surrounding the vital organs in the trunk (abdominal) area. Research shows that even if your weight and body fat remains constant, as you get older the distribution of fat changes and is more likely to shift to the trunk area especially post menopause. Ensuring you, have healthy levels of visceral fat may reduce the risk of certain diseases such as heart disease, high blood pressure, and the onset of type 2 diabetes.

The Tanita Body Composition Monitor will provide you with a visceral fat rating from 1 – 59.

**Rating from 1 – 12** 

Indicates you have a healthy level of visceral fat. Continue monitoring your rating to ensure that it stays within the healthy range.

**Rating from 13 – 59** 

Indicates you have an excess level of visceral fat. Consider making changes in your lifestyle possibly through diet changes or increasing exercise.

**Source** : Data from Columbia University (New York) & Tanita Institute (Tokyo)

**Note:**

- Even if you have a low body fat rate, you may have a high visceral fat level.
- For medical diagnosis, consult a physician.

## WHAT IS BASAL METABOLIC RATE (BMR)?

(Applicable age 18-99)

### WHAT IS BMR?

Your Basal Metabolic Rate (BMR) is the minimum level of energy your body needs when at rest to function effectively including your respiratory and circulatory organs, neural system, liver, kidneys, and other organs. You burn calories when sleeping.

About 70% of calories consumed every day is used for your basal metabolism. In addition, energy is used when doing any kind of activity however; the more vigorous the activity is the more calories are burned. This is because skeletal muscle (which accounts for approximately 40% of your body weight) acts as your metabolic engine and uses a large amount of energy. Your basal metabolism is greatly affected by the quantity of muscles you have, therefore increasing your muscle mass will help increase your basal metabolism.

By studying healthy individuals, scientists have found that as people age, their metabolic rate changes. Basal metabolism rises as a child matures. After a peak at the age of 16 or 17, it typically starts to decrease gradually.

Having a higher basal metabolism will increase the number of calories used and help to decrease the amount of body fat. A low basal metabolic rate will make it harder to lose body fat and overall weight.

### HOW DOES A TANITA BODY COMPOSITION MONITOR CALCULATE BMR?

The basic way of calculating Basal Metabolic Rate BMR is a standard equation using weight and age. Tanita has conducted in-depth research into the relationship of BMR and body composition giving a much more accurate and personalized reading for the user based on the impedance measurement. This method has been medically validated using indirect calorimetry (measuring the breath composition).\*

\*Reliability on equation for Basal Metabolic Rate:

At 2002 Nutrition Week: A Scientific and Clinical Forum and Exposition

Title: International Comparison: Resting Energy Expenditure Prediction Models:

The American Journal of Clinical Nutrition.

## WHAT IS DAILY CALORIE INTAKE (DCI)?

(Applicable age 18-99)

"Daily Calorie Intake (DCI)" is the sum of calories for basal metabolism, daily activity metabolism (activities including daily household chores), and diet-induced thermogenesis (energy used in connection with digestion, absorption, metabolism, and other eating activities). It is an estimate of how many calories you can consume within the next 24 hours to maintain your current weight.

### HOW DOES A TANITA BODY COMPOSITION MONITOR CALCULATE DCI?

$DCI = BMR \times \text{Activity Level}$

Activity Level

	1	2	3
Female	1.56	1.64	1.82
Male	1.55	1.78	2.10

Source: World Health Organization (WHO)

## WHAT IS METABOLIC AGE?

(Applicable age 18-99)

**This feature calculates your BMR and indicates the average age associated with the type of metabolism.**

If your BMR Age is higher than your actual age, it is an indication that you need to improve your metabolic rate. Increased exercise will build healthy muscle tissue, which will improve your metabolic age.

You will obtain a reading between 12 and 90. Under 12 will be displayed as "12" and over 90 displayed as "90".

# WHAT IS MUSCLE MASS?

(Applicable age 18-99)

This feature indicates the weight of muscle in your body. The muscle mass displayed includes the skeletal muscles, smooth muscles (such as cardiac and digestive muscles) and the water contained in these muscles.

Muscles play an important role as they act as an engine in consuming energy. As your muscle mass increases, your energy consumption increases helping you reduce excess body fat levels and lose weight in a healthy way.

## WHAT IS PHYSIQUE RATING?


This feature assesses your physique according to the ratio of body fat and muscle mass in your body.

As you become more active and reduce the amount of body fat, your physique rating will also change accordingly. Even though your weight may not change, your muscle mass and body fat levels may be changing making you healthier and at lower risk of certain diseases.

Each person should set their own goal of which physique they would like and follow a diet and fitness programme to meet that goal.

Result	Physique Rating	Explanation
1	Hidden Obese	<b>Small Frame Obese</b> This person seems to have a healthy body type based on physical appearance; however, they have a high body fat % with low muscle mass level.
		<b>Medium Frame Obese</b> This person has a high body fat percentage, with a moderate muscle mass level.
2	Obese	<b>Large Frame Obese</b> This person has both a high body fat % and a high muscle mass.
		<b>Low Muscle &amp; Average Body Fat %</b> This person has an average body fat % and a less than average muscle mass level.
3	Solidly-built	<b>Ave. Muscle &amp; Ave. Body Fat %</b> This person has average levels of both body fat and muscle mass.
		<b>High Muscle &amp; Ave. Body Fat % (Athlete)</b> This person has an average body fat % and higher muscle mass level than the average.
4	Under exercised	<b>Low Muscle &amp; Low Fat</b> Both body fat % and muscle mass are lower than the average.
		<b>Thin and muscular (Athlete)</b> This person has lower normal body fat % than average while having adequate muscle mass.
5	Standard	<b>Very Muscular (Athlete)</b> This person has lower normal body fat % than average, while having above-average muscle mass.
6	Standard Muscular	
7	Thin	
8	Thin and muscular	
9	Very Muscular	

Source: Data from Columbia University (New York) & Tanita Institute (Tokyo)

## WHAT IS BONE MASS?

(Applicable age 18-99)

### This feature indicates the amount of bone (bone mineral level, calcium or other minerals) in the body.

Research has shown that exercise and the development of muscle tissue are related to stronger, healthier bones. While bone structure is unlikely to make noticeable changes in a short period, it is important that you develop and maintain healthy bones by having a balanced diet and plenty of exercise. People worried about bone disease should consult their physician. People who suffer from osteoporosis or low bone densities due to advanced age, young age, pregnancy, hormonal treatment or other causes, may not get accurate estimations of their bone mass.

Below is the result of estimated bone masses of persons aged 20 to 40, who are said to have the largest amounts of bone masses, by weight. (Source: Tanita Body Weight Science Institute)

Please use the below charts as a guide to compare your bone mass reading.

#### Women: Average of estimated bone mass

Weight (lb)		
Less than 110 lb	110 lb - 165 lb	165 lb and up
4.3 lb	5.3 lb	6.5 lb

Weight (kg)		
Less than 50 kg	50 kg - 75 kg	75 kg and up
1,95 kg	2,40 kg	2,95 kg

#### Men: Average of estimated bone mass

Weight (lb)		
Less than 143 lb	143 lb - 209 lb	209 lb and up
5.9 lb	7.3 lb	8.1 lb

Weight (kg)		
Less than 65 kg	65 kg - 95 kg	95 kg and up
2,66 kg	3,29 kg	3,69 kg

#### Note:

- Persons described below may obtain varying readings and should take the values given for reference purposes only.
  - Elderly persons
  - Women during or after menopause
  - People receiving hormone therapy
- "Estimated bone mass" is a value estimated statistically based on its correlation with the fat-free amount (tissues other than the fat). "Estimated bone mass" does not give a direct judgment on the hardness or strength of the bones or the risks of bone fractures. If you have concerns over your bones, you are recommended to consult a specialist physician.

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**Tanita Health Equipment H.K. Ltd.**
**Manufacturer****Tanita Corporation**
**TANITA India Private Limited**